

World Leader in Cryosurgical Equipment Since 1966



#1 Choice of Dermatologists!
85% of dermatologists surveyed use Brymill's CRY-AC®
Largest selection of LN₂ cryosurgical spray tips, probes, & accessories

Tel: (800) 777-2796 or (860) 875-2460

Fax: (860) 872 2371

email: brymill@brymill.com

web: www.brymill.com

The Brymill Cry-Ac®Tracker®



The Cry-Ac[®]Tracker[®] is a device that measures the temperature of the skin when spraying Liquid Nitrogen via the Cry-Ac[®] delivery system to the selected skin lesion using the latest Infrared Sensoring technology.

The operator can set the Cry-Ac®Tracker® to indicate when a predetermined freeze temperature has been achieved at the lesion.

For More Precise & Effective Cryosurgery

- Consistently achieve desired skin temperature & freeze time.
- More accurately perform cryosurgery.
- Know how long to spray & when to stop.
- Quickly train others to reliably do cryosurgery.
- Achieve optimum results for patients.

Temperature Monitoring

- Infrared light sensor continuously & safely monitors skin temperature at freeze site.
- Color-coded lights shine on lesion to indicate how fast skin temperature is decreasing and when desired skin temperature is reached.

The Cry-Ac[®]Tracker[®] is not designed to be used in conjunction with contact probes.

Warranty

The Tracker® section of the Cry-Ac®Tracker® is covered by a 6 month limited warranty against manufacturing defects in materials and workmanship. The warranty does not cover damage resulting from abuse, neglect improper handling or improper use.

The Cry-Ac[®] section of the Cry-Ac[®]Tracker[®] is covered by the same 3-year warranty defined in the Cry-Ac[®] Instructions For Use booklet supplied with each unit.

Features

- Measurement accuracy better then ±3 °C.
- Minimum temperature reading -70 °C.
- Sampling rate greater than 2 Hz
- Sample area at focus less than 2 mm
- 2 x AA batteries.
- LCD screen to display ACTUAL TEMPERATURE, Tmin, and TIME AT Tmin.
- POWER- MODE buttons to turn ON/OFF and select operating MODE.
- LED to identify location of infrared temperature measurement
 - o LED color changes as a function of temperature
 - Temp 30 0 °C BLUE
 - Temp 0 Tmin °C GREEN
 - Temp Tmin °C RED
- Bleep sound to count and indicate when Tmin is achieved.
- The Cry-Ac[®] section features are the same as the Cry-Ac-3[®]

Each Cry-Ac[®]Tracker[®] is supplied with

- 1 x B and 2 x C Apertures, 1 x ½" long 20g Needle Spray
- Quick Start DVD
- Carry Case
- Instructions for Use

Model B-900

The Brymill Cry-Ac® Hand Held Liquid Nitrogen Delivery System

Brymill is the world's leading manufacturer of handheld, liquid nitrogen cryosurgery devices used to treat a wide range of common skin lesions.

Our 2 units, available in 500ml and 300ml capacity, provide the ultimate in controlled application of Liquid Nitrogen to the lesion.

The No.1 choice of dermatologists, Brymill offers unsurpassed safety, versatility, and control for precise and effective cryosurgery.



New conical apertures give a more focused spray for greater accuracy.



- Slender, fingertip trigger for maximum control and ease of viewing – also rotates easily for left-hand
- Patented Safety Autovent maximizes safety by allowing internal pressure to gradually vent as the cap is unscrewed.
- Streamlined relief valve maintains a constant operating pressure to ensure consistent and accurate freezing.
- Polypropylene Cover and Collar insulates the user's hand
- Stainless steel and brass construction for long life Sturdy Delrin base for extra stability
- Sturdy Polypropylene base for extra stability

CRY-AC[®] CRY-AC-3[®]

A set of 5 spray tips A, B, 2C's & D and a 20 gauge bent spray are included with each CRY-AC® and CRY-AC®-3 for more accurate, controlled treatment of lesions of different sizes and types.

Model	Description	Height	Weight		Static Holding Time
			Full	Empty	
B-700	CRY-AC® (500ml)	11.25in	31.6oz	19.0oz	20 – 24 hours
		(286mm)	(895g)	(540g)	
B-800	CRY-AC®-3 (300ml)	8.5in	24oz	14.5oz	10 – 12 hours
	, ,	(215mm)	(618g)	(412g)	

Starter Packs

Our starter pack includes

- * Cry-Ac®-3 (or Cry-Ac®) unit with 5 spray tips and a 20g Bent spray (See above)
- * A 20 liter Dewar See details below
- * Method of Liquid Withdrawal see details below



Image of BRY-1001

Family Practice Package

MODEL	DESCRIPTION
	Cry-Ac®-3 – 10oz., Withdrawal Tube, 20SC Dewar, Cryoplate (see page 5)
BRY-1001	Cry-Ac® – 16oz., Withdrawal Tube, 20SC Dewar, Cryoplate (see page 5)

Dermatology Package

MODEL	DESCRIPTION
BRY-1002	Cry-Ac®-3 – 10oz., Withdrawal Device, 20 liter Dewar, Cryoplate (see page 5)
BRY-1003	Cry-Ac® – 16oz., Withdrawal Device, 20 liter Dewar, Cryoplate (see page 5)

Liquid Nitrogen Storage & Withdrawal Equipment



		Static Weight Weight		Weight	Dimer	nsions
Size	Model	Holding	Empty	Full	Overall	Outside
		Time			Height	Diameter
5 Litre	501-5	4 – 5 wks.	8 lbs.	17 lbs.	18.2in.	8.8in.
0 2 0	55.5	· · · · · · · · · · · · ·	0 .50.		462mm	222mm
10 Litre	501-10	6 – 8 wks.	13 lbs.	30 lbs.	21.5in.	10.3in.
		o o men		00 1001	546mm	260mm
20 Litre	501-20	8 – 12 wks.	19 lbs.	53 lbs.	24.7in.	14.5in.
					627mm	368mm
20 Litre	501-20SC	Long Last	19 lbs.	53 lbs.	24.7in	14.5in
		200 days			627mm	368mm
30 Litre	501-30	14 – 16 wks.	26 lbs.	77 lbs.	24.1in.	17in.
2.0.0					611mm	432mm
50 Litre	501-50	14 – 17 wks.	34 lbs.	123 lbs.	30.7in.	17in.
JULIA DE LINO			J . 1651	.20 100.	79mm	432mm



Model
Description
Withdrawal Device (*) - For Dewars of 20 litres and over

Withdrawal Tube (*) - For Dewars of 5 up to 20 litres

Replacement Filter for Withdrawal Device

Replacement Filter for Withdrawal Tube



Withdrawal Device



Roller Base (*)
Recommended for Dewars
of 20 ltrs and over



Dewar Cap (*)
Replacement for worn
or damaged Dewar
Caps

Note: When ordering the asterisked items, please specify Dewar manufacturer, model # and capacity

Model	Description	Mid-arm length: 14 - 15in
605-S	Cry-Gloves® - Small	The most popular length of cryogenic glove offers maximum thermal protection for your
605-M	Cry-Gloves® - Medium	hands and forearms when working in hazardous, ultra-cold environments.
605-L	Cry-Gloves® - Large	CAUTION Not intended for immersion in liquid nitrogen or other cryogenic liquids.

Open Sprays




106-20	
107-16	
107-18	
107-20	
Model	Description
106	Original Straight Spray (not pictured)
106-20	3in. x 20g Straight Spray
106-18	3in. x 18g Straight Spray (not pictured)
106-16	3in. x 16g Straight Spray (not pictured)
107-16	1.5in. x 16g Straight Spray Needle
107-18	1.5in. x 18g Straight Spray Needle
107-20	1.5in. x 20g Straight Spray Needle

Spray Extensions & Adapters



301	308	309	310-45 & 310
Model	Description		
301	Back Vent Adapter		
308	Luer Lock Adapter		
309	Malleable Extension (3, 6, 9 & 12")		
310	Right Angle Adapter (Recommended for use with Cryo Chambers)		

Cryoplate & Cryocones



303	Cryoplate		
A transparent lex			
various sized co	various sized conical openings (3, 5,		
8 and 10mm). Used with apertures,			
the cryoplate provides even greater			
localization of freezing and			
protection of sensitive areas such as			

Description



Model	Description	
304	Cryocones	
A set of 6 varied size		
cones used to conce		
within a limited area. Especially		
helpful for irregular shaped lesions		
because they can be	e shaped to the	
lesion.		
Sizes are 6, 11, 16,	25, 30 and	
38mm.		

Closed Probes

Mini Probes



Model	Description
214-1	1mm Mini Probe
214-2	2mm Mini Probe
214-3	3mm Mini Probe
214-4	4mm Mini Probe
214-5	5mm Mini Probe
214-6	6mm Mini Probe
214-SP	Sharp Pointed Mini Probe
214-15	1.5x3mm Rectangular Mini Probe
214-25	2.5x5mm Rectangular Mini Probe

Conical Probes



Model	Description
216	Sharp-pointed Conical Probe
203-1	1mm Conical Probe
203-2	2mm Conical Probe
203-3	3mm Conical Probe
203-4	4mm Conical Probe
203-5	5mm Conical Probe
203-6	6mm Conical Probe

Ball Probes



Model	Description
201-6	6mm Ball Probe
201-8	8mm Ball Probe
201-1	1cm Ball Probe
201-2	2cm Ball Probe
201-3	3cm Ball Probe

Flat Probes



Model	Description
205-1	1cm Flat Probe
205-2	2cm Flat Probe
205-2.5	2.5cm Flat Probe
205-3	3cm Flat Probe
205-8	8mm Gentle Probe

Cervical Probes



Model	Description
206	19mm Domed Cervical Probe
207	19mm Flat Cervical Probe
208	25mm Domed Cervical Probe
209	25mm Flat Cervical Probe
210	19x5x5mm Cervical Probe
212	25x5x5mm Cervical Probe

Miscellaneous Probes


Model	Description
202	1cm Concave Probe
213	6 x 32mm Advanced Hemorrhoid Probe

Additional Products



Model	Description
601	Cryotote

Durable insulated bag with long handle for convenient hands-free carrying of CRY-AC® or CRY-AC®-3 off-site



Cutaneous Cryosurgery: Principles & Clinical Practices: - Third Edition.

The leading textbook that covers essential information for starting and maintaining a cryosurgical practice. Describes diagnoses and cryosurgical treatments for many benign, pre-malignant, and malignant lesions. Extensively illustrated with numerous color photos that can aid in diagnosis and treatment of various dermatological conditions. By Drs. Rodney Dawber, Graham Colver, and Arthur Jackson.

Available from your Medical Book Supplier

ISBN number: - 1084184-552-3



	cription
350 Cryo	Tweezers

The Brymill Cryo Tweezers have been designed for the localized treatment of Skin Tags. By using the Brymill Cryo Tweezers you can control freezing easily and accurately.

Supplied with a Delrin Cup Holder and 10 Foam Cups.

To reorder more Styrofoam Cups simply Google Dart 4J4



Model	Description
360	Histology Package

The Brymill Hystology Package consists of a 10 inch Malleable Spray Extension and a Wide Spray (Acne) Aperture. Commonly used in MOHS Labs.

Quality Assurance Guaranteed

All products designed and manufactured by Brymill Cryogenic Systems are done so in accordance with appropriate F.D.A, Canada Health and European Regulatory requirements, underscoring Brymill's commitment to supplying safe, quality products.









Brymill Offices



Brymill Cryogenic Systems

World Headquarters

105 Windermere Avenue Ellington, CT 06029 Tel: (860) 875 2460

Email: brymill@brymill.com Web: www.brymill.com

Southern California Office

Tel: (858) 586 0058 Fax: (858) 586 0081 Email: brymill@brymill.com Web: www.brymill.com



Brymill Cryogenic Systems (UK) Ltd

54 Mayfield Ridge, Hatch Warren
Basingstoke. RG22 4RS
Tel: +44 (0)1256 841045 Fax: +44 (0)1256 350330

Email: sales@brymilluk.com Web: www.brymill.co.uk



Brymill Cryogenic Systems – Latin America

G-Mistral 3174
Buenos Aires
1419, Argentina
Tel/Fax: 54-11-45714557
Email: crusso@brymill.com
Web: www,brymill.com/espanol

Brymill's US Dealers and International Distributors

Fax: (860) 872 2371

To find your nearest US Dealer, International Distributor or Repair Centre visit our web site at http://www.brymill.com/distributors repair centres.html

Why Brymill

- World Leader Since 1966
- #1 Choice of 85% of Dermatologist
- Largest Selection of LN2 Equipment, Tips Probes & Accessories
- Superior Quality, Reliability & Durability
- 3-Year Limited warranty against manufacturing defects
- 30-Day Unconditional Money-back Guarantee

What Dermatologists Say...

"Brymill hand-held units are an indispensable part of my practice. Because of the sizes and design, they are comfortable to hold and easy to use. They function excellently, without clogging, and seem to last forever. The large selection of accessories available makes it possible for me to treat many types of benign, pre-malignant and malignant lesions. I would recommend these units to any physician the performing cryosurgery."

Emanuel G. Kuflik, M.D., Adult & Child Dermatology Center, Lakewood, NJ

Author of Cryosurgical Treatment for Skin Cancer

"CRY-AC® is what dermatologists use every day in their practice.... I especially like Brymill's unit because it is very sturdy and holds up well after years of heavy use. They stand behind their units. Brymill is a company you can trust."

Gloria F. Graham, M.D., Grahams' Dermatology Services, P.A., Winston-Salem, NC Author of cryosurgery chapter in R. Wheeland's Cutaneous Surgery

What other Physicians Say...

"Our family practice group has been using Brymill's CRY-AC® for about 15 years. We have achieved superb results with our patients. Brymill has constantly provided us with excellent service and a very personalized touch. I recommend the Brymill system without hesitation."

Neil H. Brooks, M.D., Past President, AAFP, Health Wise Medical Associates, Vernon, CT

"I use my Brymill CRY-AC® 20-30 times per week for treating everything from warts to molluscum to AK's to skin tags. It's incredibly easy to use. I trained my nurse practitioner to use it, and she uses it as often as I do." Mark Moores, M.D., Webster Family Health Center, Dudley, MA

"Brymill's CRY-AC® is a very physician-friendly device for a busy family practitioner's office. All five doctors in our office use it, and at least one of the five of us uses it daily. The CRY-AC® quickly pays for itself many times over."

H. Andrew Selinger, M.D., Pro Health Physicians Family Medical Group, Bristol, CT