

EASYPULSE TOC

Total Oxygen Concentrator







The Perfect Balance of Portable and Stationary Oxygen Delivery





Total Oxygen Concentrator

Model: PM4400 Series

Total Oxygen Concentrator

It is the perfect balance of mobile and stationary oxygen delivery. It delivers continuous flow and pulse flow all in the same unit weighing less than 11.5 pounds. No need to have 2 separate oxygen concentrators.

Easy to Use

The EasyPulseTOC has an easy to use patient interface. One simple knob easily sets the oxygen flow to either pulse or continuous mode. Unlike many other TOC's, there are no confusing electronic controls.

Perfect for the Non-Delivery Oxygen Model

The EasyPulseTOC is the ideal product to help you reduce costly cylinder deliveries and eliminate the need for stationary oxygen concentrators.



Makes Mobility Easy

The convenient fold-up cart, with its adjustable telescoping handle and oversized wheels, is easily attached to the TOC. This enables the patient to take the unit anywhere and by any means. It is FAA approved for use on airlines and fits easily under an airline seat.

All You Need is in the Box

Open the box and you are set to go! It contains the Concentrator, Cart, AC Power Supply, AC Power Cord, DC Power Cord, Battery, Accessory Bag and User Manual.

Excellent Battery Life

The battery lasts over 3 hours. And with optional, easy to change extra batteries, you can extend your mobility as long as you wish.

Specifications

Dimensions Height 14.63 in. (37.16 cm)

Width 10.50 in. (26.67 cm) Depth 7.00 in. (17.78 cm)

Weight

11.5 lbs (5.2 kg)

Operating Temperature

TOC: 41° to 104°F (5° to 40°C)

Charging Battery: 32° to 113°F (0° to 45°C)

Altitude

Sea Level to 10,000 ft (3048 m)

Oxygen Concentration

87 to 95% at all settings

Noise Level

50.6 dBa on Setting 2 Continuous 42.9 dBa on Pulse Setting 1

Battery Duration

| Pulse Mode | Continuous Mode | |
|----------------------|-----------------|--|
| Setting 1: 3.3 Hours | 1/8 5.0 hours | |
| Setting 2: 3.0 Hours | 1/4 4.3 hours | |
| Setting 3: 2.8 Hours | 1/2 3.3 hours | |
| Setting 4: 2.5 Hours | 1 1.8 hours | |
| Setting 5: 2.3 Hours | 2 0.8 hours | |

Cannula Requirement

Pulse Mode: Maximum 7 foot long adult

standard single lumen

oxygen cannula

Continuous Mode: Maximum 30 foot long adult

standard single lumen

oxygen cannula

Specifications are subject to change without prior notice.

Ordering Information

| Part Number | Description |
|-------------|--|
| PM4400 | EasyPulse TOC Total Oxygen Concentrator - includes AC and DC Power Supplies, Cords, Carry Bag, Cannula & Battery |
| PM4400EB | EasyPulse TOC Total Oxygen Concentrator - includes AC and DC Power Supplies, Cords, Carry Bag, Cannula & 2 Batteries |

In the Box



Desktop Charger

Model: 508306



Batteries sold separately

Specifications

Charger Base

Dimensions

Height 3.00 in. (7.62 cm) Width 6.20 in. (15.75 cm) Depth 7.80 in. (19.81 cm

Weight

0.80 lbs (0.36 kg)

AC Power Supply

Weight

0.77 lbs (0.35 kg)

Input

100-240V ~ 50-60HZ 1.6A

Specifications are subject to change without prior notice.

Charge Time

Charging a <u>Single</u> fully depleted Battery will take approximately 2 hrs. Charging <u>Two</u> fully depleted Batteries will take approximately 4 hrs.

Accessories

Desktop Charger



Batteries sold separately

Ordering Information: 508306 **Desktop Charger**

Disposable Humidifier



Ordering Information:

PM500 Disposable Humidifier (case of 50)

Reusable Humidifier



Ordering Information:

PM510 Reusable Humidifier

Air Inlet Filter



Ordering Information: 508186 Air Inlet Filter

Battery



Ordering Information:

508016 **Battery**

Humidifier Connector with Tube



Ordering Information:

Humidifier Connector 508333 with Tube

Transportable Cart



Ordering Information:

508260 Transportable Cart

Accessory Bag



Ordering Information:

508301 Accessory Bag

Desktop Charger Base



Ordering Information:

508367 Desktop Charger Base

AC Power Supply





Ordering Information:

506629

508188 **AC Power Supply** (PM4400)

> **AC Power Supply** (Desktop Charger)

AC Power Cord





Ordering Information:

AC Power Cord 508291 (PM4400) **AC Power Cord** 506961 (Desktop Charger)

DC Power Cord



Ordering Information:

DC Power Cord 508189 (PM4400 only)

PRECISION MEDICAL.

Precision Medical, Inc.

300 Held Drive • Northampton, PA 18067 USA

Phone: 610-262-6090 Toll-Free Phone: 800-272-7285 Fax: 610-262-6080 Toll-Free Fax: 800-352-1240

www.precisionmedical.com